

ICE WRAP INSTRUCTIONS

NEOPRENE WRAP WITH ICE BAG

- Insert Ice Bag through hole (Screw cap opening should be same side as SGA Label)
- Ice Bag - Unscrew lid insert ice to require amount (about 1/2 to 2/3 is ideal)
- Screw Lid back on firm but not too tight
- Place Ice Bag on swollen area then hold in place with Neoprene Wrap using Velcro to reposition for comfort

IMPORTANT:

Never Ice for more than 15-20 minutes to avoid frost bite

Before re icing swollen area ensure it has no numbness (Suggested is no sooner than 1 hour)

REMEMBER:

Rest, Compression and elevation assist the healing process and minimise damage to the affected area

To view more information on our products or to Buy Online Visit us at WWW.SUPPORTSGUARDS.COM

SUPPORT GUARDS AUSTRALIA